Thermofluide / Heat Transfer Fluids

Inspired by **temperature**

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Huber heat transfer fluids have the best thermodynamic and environmentally friendly characteristics. The choice of the correct fluid is decisive for efficient heat transfer. The working temperature range is generally the most important criterion for choosing a particular fluid. The use of the correct heat transfer fluid guarantees reliable and safe operation, optimum temperature control and a long life. The heat transfer fluids in this brochure have been chosen for use with Huber temperature control units.

If you require help in choosing the heat transfer fluid, we are happy to help:

Hotline Sales: +49 781 9603 123

Hotline Technical Support: +49 781 9603 244

Safety data sheets and technical information for our heat transfer fluids can be found under www.huber-online.com

Tip for Unistats:

For Unistats, we recommend the heat transfer fluid "DW-Therm" which has a high specific heat capacity and so contributes towards an optimum heat transfer. In addition, DW-Therm covers a large temperature range, and has an impressively long service life.

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Drain tap for easy filling with thermofluid directly from the canister.

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
The table shows an indicative overview of the heat transfer fluids that may be used in each unit. When selecting the thermal fluids, the operating temperatures, the operating instructions and the data sheet of the temperature control unit as well as application-specific features must be observed.

DW-Therm M90.200.02
DW-Therm HT P20.340.32
SilOil P20.225.275.50
SilOil M20.195.235.20
SilOil M40.165.220.10
SilOil M60.115.200.05
SilOil M80.055.03
SilOil M80.100.250.03
SilOil M90.055.170.02
SynOil M10.120.08
MinOil P20.190.40
MEG
Wasser/Water

Unistat Temperiersysteme / Temperature Control Systems	
Unistat Petite Fleur, Grande Fleur, Tango – 430w	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Unistat 510 – 540w	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Unistat 610 – 640w	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Unistat 645 – 680w	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Unistat 705 – 825w	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Unistat 904 – 950w	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Unistat 1005 – 1015w	auf Anfrage / on request
Unistat T305 – T402, TR401 – TR402, Chili	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Unimotive	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Umwälzkühler / Chillers	
Piccolo	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Minichiller	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Unichiller 015 – 025	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Unichiller P007 – P025	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Unichiller 017T – 500T	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Unichiller 050 – 230, P050 – P100w	● ● ● ● ● ● ● ● ● ● ● ● ● ●
RotaCool	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Eintauchkühler / Immersion Cooler TC45 – TC100	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Bad- und Umwälzthermostate / Bath Circulators	
Einhänge-Thermostate / Immersion Circulators	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Badthermostate, Polycarbonat / Bath Circulators, Polycarbonate	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Bad-/Umwälzthermostate, Edelstahl / Bath Circulators, Stainless Steel	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Visco-Thermostate / Visco Baths	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Brücken-Thermostate / Bridge Circulators	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Kältethermostate / Cooling Circulators	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Ministat	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Variostat	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Specials	
Bier-Forciertest-Thermostat / Beer Force-Ageing-Test-Bath	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Hotbox	● ● ● ● ● ● ● ● ● ● ● ● ● ●
Heat Transfer Station	● ● ● ● ● ● ● ● ● ● ● ● ● ●

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


DW-Therm is a mixture of isomeric triethoxysilanes and has been developed for hydraulically sealed systems.

We recommend the heat transfer fluid DW-Therm for Unistats. DW-Therm has a high specific heat capacity and so contributes towards an optimum heat transfer. In addition, DW-Therm covers a large temperature range, and has an impressively long service life.

Eigenschaften <i>Properties</i>	Wert <i>Value</i>
Arbeitstemperatur °C <i>Working temperature °C</i>	-90 ... 200
Flammpunkt °C <i>Flash Point °C</i>	101
Brennpunkt °C <i>Fire Point °C</i>	112
Viskosität mm ² /s (kinematisch bei 25 °C) <i>Viscosity mm²/s (kinematic at 25 °C)</i>	2
Dichte g/cm ³ (bei 20 °C) <i>Density g/cm³ (at 20 °C)</i>	0,879
Stockpunkt °C <i>Pour Point °C</i>	-137
Zündtemperatur °C <i>Ignition temperature °C</i>	265
Farbe <i>Colour</i>	leicht gelblich, klar / <i>yellowish, clear</i>
Thermischer Ausdehnungskoeffizient 10 ⁻⁵ /K <i>Thermal expansion coefficient 10⁻⁵/K</i>	96,98
Wärmeleitfähigkeit W/(m·K) bei 30 °C <i>Heat conductivity W/(m·K) at 30 °C</i>	0,116
Materialunverträglichkeit <i>Material incompatibility</i>	Kupfer, Aluminium, Buntmetalle / <i>copper, aluminum, non-ferrous metals</i>

Please pay attention to our safety data sheets and technical data sheets before use. These documents can be downloaded from www.huber-online.com or requested directly from us.

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DW-Therm HT is a mixture of partially hydrogenated naphthalenes. It is especially suited for high temperature applications using hydraulically sealed thermostats.

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
SilOil P20.225/275.50 is a low viscosity silicon oil, that due to its special characteristics is preferred as a cooling as well as a heating medium which can be used in temperature control units as well as in heating systems.

SilOil P20.225/275.50 may be used within the range 20 °C to 225/275 °C in connection with externally closed systems (Unistats). It can be used in open bath systems, however should be noted that P20.225/275.50 when used at high temperatures, in the presence of an oxidising medium such as air, or substances which act as a catalyst such as acids, alkalis, and various metallic compounds, can alter chemically. In the presence of oxidation products, an increase in viscosity or even "gelling" of the oil through a polymerisation reaction is to be expected, while in contact with catalysing substances mostly a depolymerisation process and a resulting lowering of viscosity can be observed.

Eigenschaften <i>Properties</i>	Wert <i>Value</i>
Arbeitstemperatur °C <i>Working temperature °C</i>	20 ... 225/275
Flammpunkt °C <i>Flash Point °C</i>	>300
Brennpunkt °C <i>Fire Point °C</i>	keine Angabe / <i>not specified</i>
Viskosität mm ² /s (kinematisch bei 25 °C) <i>Viscosity mm²/s (kinematic at 25 °C)</i>	50
Dichte g/cm ³ (bei 20 °C) <i>Density g/cm³ (at 20 °C)</i>	0,96
Stockpunkt °C <i>Pour Point °C</i>	-55
Zündtemperatur °C <i>Ignition temperature °C</i>	450
Farbe <i>Colour</i>	farblos, klar / <i>colourless, clear</i>
Thermischer Ausdehnungskoeffizient 10 ⁻⁵ /K <i>Thermal expansion coefficient 10⁻⁵/K</i>	95
Wärmeleitfähigkeit W/(m·K) bei 50 °C <i>Heat conductivity W/(m·K) at 50 °C</i>	0,14

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


P20.225/275.50 is chemically inert. Contact with concentrated alkaline solutions, oxidising acids, chlorine gas and hydrochloric acid should be avoided, particularly at elevated temperatures. Incompatible with silicone rubber.

Please pay attention to our safety data sheets and technical data sheets before use. These documents can be downloaded from www.huber-online.com or requested directly from us.

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
SilOil M20.195/235.20 is a low-viscosity silicone fluid which, as a result of its special property profile, is particularly suitable for use as a cold and heat transfer medium in cryostats, thermostats and heat transfer installations.

SilOil M20.195/235.20 can be used in the range from -20 °C to 195 °C (for open systems) and to 235 °C in externally sealed systems (Unistats). It should be borne in mind that, at high temperatures, SilOil M20.195/235.20 can be chemically altered by oxidising media, such as air, or substances with a catalytic effect, such as acids, lyes and various metal compounds. An increase in viscosity, and possibly even gelling of the fluid owing to crosslinking reactions, must be expected in the presence of oxidising agents, while contact with products having a catalytic effect usually induces a process of depolymerisation, resulting in a drop in viscosity.


Eigenschaften <i>Properties</i>	Wert <i>Value</i>
Arbeitstemperatur °C <i>Working temperature °C</i>	-20 ... 195/235
Flammpunkt °C <i>Flash Point °C</i>	> 200
Brennpunkt °C <i>Fire Point °C</i>	> 290
Viskosität mm ² /s (kinematisch bei 25 °C) <i>Viscosity mm²/s (kinematic at 25 °C)</i>	20
Dichte g/cm ³ (bei 20 °C) <i>Density g/cm³ (at 20 °C)</i>	0,95
Stockpunkt °C <i>Pour Point °C</i>	< -70
Zündtemperatur °C <i>Ignition temperature °C</i>	> 400
Farbe <i>Colour</i>	farblos, klar / <i>colourless, clear</i>
Thermischer Ausdehnungskoeffizient 10 ⁻⁵ /K <i>Thermal expansion coefficient 10-5/K</i>	101
Wärmeleitfähigkeit W/(m·K) bei 25 °C <i>Heat conductivity W/(m·K) at 25 °C</i>	0,151
Materialunverträglichkeit <i>Material incompatibility</i>	Silikonkautschuk / <i>Silicone rubber</i>

Please pay attention to our safety data sheets and technical data sheets before use. These documents can be downloaded from www.huber-online.com or requested directly from us.

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
SilOil M40.165/220.10 is a low-viscosity silicone fluid which, as a result of its special property profile, is particularly suitable for use as a cold and heat transfer medium in cryostats, thermostats and heat transfer installations.


SilOil M40.165/220.10 can be used in the range from -40 °C to 165 °C and to 220 °C in externally sealed systems (Unistats).

It should be borne in mind that, at high temperatures, SilOil M40.165/220.10 can be chemically altered by oxidising media, such as air, or substances with a catalytic effect, such as acids, lyes and various metal compounds. An increase in viscosity, and possibly even gelling of the fluid owing to crosslinking reactions, must be expected in the presence of oxidising agents, while contact with products having a catalytic effect usually induces a process of depolymerisation, resulting in a drop in viscosity.


Eigenschaften <i>Properties</i>	Wert <i>Value</i>
Arbeitstemperatur °C <i>Working temperature °C</i>	-40 ... 165/220
Flammpunkt °C <i>Flash Point °C</i>	> 170
Brennpunkt °C <i>Fire Point °C</i>	> 230
Viskosität mm ² /s (kinematisch bei 25 °C) <i>Viscosity mm²/s (kinematic at 25 °C)</i>	10
Dichte g/cm ³ (bei 20 °C) <i>Density g/cm³ (at 20 °C)</i>	1
Stockpunkt °C <i>Pour Point °C</i>	< -90
Zündtemperatur °C <i>Ignition temperature °C</i>	> 400
Farbe <i>Colour</i>	farblos, klar / <i>colourless, clear</i>
Thermischer Ausdehnungskoeffizient 10 ⁻⁵ /K <i>Thermal expansion coefficient 10-5/K</i>	103
Wärmeleitfähigkeit W/(m·K) bei 25 °C <i>Heat conductivity W/(m·K) at 25 °C</i>	0,142
Materialunverträglichkeit <i>Material incompatibility</i>	Silikonkautschuk / <i>Silicone rubber</i>

Please pay attention to our safety data sheets and technical data sheets before use. These documents can be downloaded from www.huber-online.com or requested directly from us.

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

SilOil M60.115/200.05 is a low-viscosity silicone fluid which, as a result of its special property profile, is particularly suitable for use as a cold and heat transfer medium in cryostats, thermostats and heat transfer installations.

SilOil M60.115/200.05 can be used in the range from -60 °C to 115 °C (for open systems) and to 200 °C in externally sealed systems (Unistats). It should be borne in mind that, at high temperatures, SilOil M60.115/200.05 can be chemically altered by oxidising media, such as air, or substances with a catalytic effect, such as acids, lyes and various metal compounds. An increase in viscosity, and possibly even gelling of the fluid owing to crosslinking reactions, must be expected in the presence of oxidising agents, while contact with products having a catalytic effect usually induces a process of depolymerisation, resulting in a drop in viscosity.


Eigenschaften <i>Properties</i>	Wert <i>Value</i>
Arbeitstemperatur °C <i>Working temperature °C</i>	-60 ... 115/200
Flammpunkt °C <i>Flash Point °C</i>	> 120
Brennpunkt °C <i>Fire Point °C</i>	> 160
Viskosität mm ² /s (kinematisch bei 25 °C) <i>Viscosity mm²/s (kinematic at 25 °C)</i>	5
Dichte g/cm ³ (bei 20 °C) <i>Density g/cm³ (at 20 °C)</i>	0,92
Stockpunkt °C <i>Pour Point °C</i>	< -100
Zündtemperatur °C <i>Ignition temperature °C</i>	> 400
Farbe <i>Colour</i>	farblos, klar / <i>colourless, clear</i>
Thermischer Ausdehnungskoeffizient 10 ⁻⁵ /K <i>Thermal expansion coefficient 10-5/K</i>	108
Wärmeleitfähigkeit W/(m·K) bei 25 °C <i>Heat conductivity W/(m·K) at 25 °C</i>	0,133
Materialunverträglichkeit <i>Material incompatibility</i>	Silikonkautschuk / <i>Silicone rubber</i>

Please pay attention to our safety data sheets and technical data sheets before use. These documents can be downloaded from www.huber-online.com or requested directly from us.

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

SilOil M80.055.03 is a low-viscosity silicone fluid which, as a result of its special property profile, is particularly suitable for use as a cold and heat transfer medium in cryostats, thermostats and heat transfer installations.

The fluid can be used in the range from -80 °C to 55 °C. It should be borne in mind that, at high temperatures, SilOil M80.055.03 can be chemically altered by oxidising media, such as air, or substances with a catalytic effect, such as acids, lyes and various metal compounds. An increase in viscosity, and possibly even gelling of the fluid owing to crosslinking reactions, must be expected in the presence of oxidising agents, while contact with products having a catalytic effect usually induces a process of depolymerisation, resulting in a drop in viscosity.


Eigenschaften <i>Properties</i>	Wert <i>Value</i>
Arbeitstemperatur °C <i>Working temperature °C</i>	-80 ... 55
Flammpunkt °C <i>Flash Point °C</i>	> 62
Brennpunkt °C <i>Fire Point °C</i>	> 110
Viskosität mm ² /s (kinematisch bei 25 °C) <i>Viscosity mm²/s (kinematic at 25 °C)</i>	3
Dichte g/cm ³ (bei 20 °C) <i>Density g/cm³ (at 20 °C)</i>	0,9
Stockpunkt °C <i>Pour Point °C</i>	< -100
Zündtemperatur °C <i>Ignition temperature °C</i>	> 400
Farbe <i>Colour</i>	farblos, klar / <i>colourless, clear</i>
Thermischer Ausdehnungskoeffizient 10 ⁻⁵ /K <i>Thermal expansion coefficient 10⁻⁵/K</i>	111
Wärmefähigkeit W/(m·K) bei 25 °C <i>Heat conductivity W/(m·K) at 25 °C</i>	0,124
Materialunverträglichkeit <i>Material incompatibility</i>	Silikonkautschuk / <i>Silicone rubber</i>

Please pay attention to our safety data sheets and technical data sheets before use. These documents can be downloaded from www.huber-online.com or requested directly from us.

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

SilOil M80.100/250.03 is a low-viscosity silicone fluid which, as a result of its special property profile, is particularly suitable for use as a cold and heat transfer medium in cryostats, thermostats and heat transfer installations.

SilOil M80.100/250.03 can be used in the range from -80 °C to 100 °C (for open systems) and to 250 °C with externally sealed systems (Unistats). It should be borne in mind that, at high temperatures, SilOil M80.100/250.03 can be chemically altered by oxidising media, such as air, or substances with a catalytic effect, such as acids, lyes and various metal compounds. An increase in viscosity, and possibly even gelling of the fluid owing to crosslinking reactions, must be expected in the presence of oxidising agents, while contact with products having a catalytic effect usually induces a process of depolymerisation, resulting in a drop in viscosity.

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

Low viscous, linear polydimethylsiloxane particularly suitable as environmentally friendly cold- and heat transfer media.




SilOil M90.055/170.02 can be used in the range from -90 °C to 55 °C (for open systems) and to 170 °C in connection with externally sealed systems (for Unistats).

SilOil M90.055/170.02 is characterised by a unique combination of low viscosity and a relatively high flash point. In contrast to competitive products, SilOil M90.055/170.02 oil is not classified as dangerous good and hazardous substance. It allows an economic use even at low temperatures due to its low viscosity. Since the SilOil M90.055/170.02 is based on polydimethylsiloxane this material is easily miscible with other media on the basis of polydimethylsiloxane. SilOil M90.055/170.02 oil is colorless, odorless and non-toxic.

Eigenschaften <i>Properties</i>	Wert <i>Value</i>
Arbeitstemperatur °C <i>Working temperature °C</i>	-90 ... 55/170
Flammpunkt °C <i>Flash Point °C</i>	> 90
Brennpunkt °C <i>Fire Point °C</i>	106
Viskosität mm ² /s (kinematisch bei 25 °C) <i>Viscosity mm²/s (kinematic at 25 °C)</i>	2
Dichte g/cm ³ (bei 25 °C) <i>Density g/cm³ (at 25 °C)</i>	0,87
Stockpunkt °C <i>Pour Point °C</i>	< -120
Zündtemperatur °C <i>Ignition temperature °C</i>	> 400
Farbe <i>Colour</i>	farblos, klar / colourless, clear
Thermischer Ausdehnungskoeffizient 10 ⁻⁵ /K <i>Thermal expansion coefficient 10⁻⁵/K</i>	125
Wärmeleitfähigkeit W/(m·K) bei 25 °C <i>Heat conductivity W/(m·K) at 25 °C</i>	keine Angabe / not specified
Materialunverträglichkeit <i>Material incompatibility</i>	Silikonkautschuk / Silicone rubber

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SynOil is a low viscosity Poly-alpha-olefin based fluid with a strongly enhanced low temperature characteristic. SynOil has excellent thermal stability as well as a very low tendency to evaporate. The fluid has a high viscosity index, which enables the product to flow, even at low temperatures, and not to build a disruptive film at high temperatures.

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

P20.190.40 is a high quality mineral oil based heat transfer fluid. It meets the technical safety requirements and tests according to DIN 4754 as well as the recommendations of VDI 3033 (Construction, operation and maintenance of heat transfer systems).

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


Heat transfer fluid with frost and corrosion protection characteristics for use in technical heating and cooling systems. As basis for the freezing protection, Monoethylene Glycol (MEG) is used, which due to its very high boiling point minimises loss through evaporation.

The MEG is diluted according to the table on the right side with water. Pure MEG is unsuitable for temperature control applications with Huber units.

Use of water / monoethylene glycol (MEG) mixtures in a temperature control unit:

The water for the mixture should have a chloride content of < 50 mg/l.. Please refer to the table (page 30) for the specifications of the water content.

Our monoethylene glycol contains corrosion inhibitors. The minimum mixing ratio for sufficient corrosion protection is 20 % MEG. This corresponds to an application temperature of down to -10 °C. The mixing ratio should not exceed 50 % MEG (viscosity). This corresponds to a pour point temperature of approx. -37 °C.

Exceptions to this are temperature control units specially suited for higher mixing ratios, such as: The Unimotive series. With these, mixing ratios of up to 60 % MEG can be used.

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* depends on the mixing ratio and the temperature control unit

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A protective additive for prevention and fighting the build up of algae in open bath units. It is used with water and prevents the build up of algae, bacteria and micro organisms by preventing germination.

Bestellnummer Order Number	Menge Quantity
6172	0,1 Ltr.

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Water is often used as a thermal fluid or cooling water, e.g. in temperature control units. To avoid corrosion and thus damage to the unit, the water used must have certain quality characteristics/limit values.

The limit values listed in the following table must be complied with if water is to be used as thermofluid or cooling water in temperature control units*. (Before using water as a thermal fluid, please check whether the temperature control unit may be operated with water).

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This resistance table is not complete and serves only as a guide. The values stated in it are guideline values which may have to be restricted under certain conditions. Any corrosion that may occur is a very complex process and is triggered by various ingredients, often in combination.

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Our recommendation:

- Colour, odour and particle content: clear/colourless, without odour and no visible particles/turbidity.
- Calcium carbonate: 0.7 mmol/l ... 1.5 mmol/l
- For purified, distilled and pure water and ultrapure water: add 0.1 g/l soda (Na_2CO_3)

The following describes how to maintain water quality and the forms and sources of corrosion that occur in the area of temperature control equipment.

Maintaining water quality:

- If only tap water is available, the content of chloride ions must be as low as possible to avoid corrosion.
- If the water circuit is operated at temperatures below 50 °C and under the influence of light, a biocide can be used to prevent algae growth (do not use biocides containing chloride!)
- If air pollutants or organic compounds are introduced into the water circuit, change the water regularly
- Alloyed standard stainless steels such as V2A (St. 1.4301) are not resistant to corrosion caused by chloride ions.
- It is not possible to „neutralise“ the chloride ions
- If an acidic chloride-containing solution is neutralised, e.g. by adding sodium carbonate solution, nothing changes in the concentration of the chloride ions. The corrosive effect remains

Differentiation of water qualities in terms of designation and conductivity

The quality of water is usually subdivided according to electrical conductivity (in $\mu\text{S}/\text{cm}$). There are different standards and definitions of purified water, which differ slightly depending on the area of application. The following table illustrates different qualities of water.

Bezeichnung / Name	Leitfähigkeit / Conductivity [$\mu\text{S}/\text{cm}$]
Trinkwasser (weich ... hart) / Drinking water (soft ... hard)	150 ... 800
Gereinigtes Wasser, salzarmes Wasser / Purified water, low-salt water	< 50
Destilliertes Wasser / Distilled water	< 5
Reines Wasser / Pure water	< 1
Reinstwasser, hochreines Wasser / Ultra-pure water, high-purity water	< 0,1

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The main cause of corrosion: chloride ions

Chloride ions are the main cause of corrosion when water is used as a thermal fluid in contact with metallic components. They are particularly responsible for critical types of corrosion such as pitting and stress corrosion cracking. Chloride ions are adsorbed by (steel) surfaces and destroy or prevent the formation of a protective (passivating) oxide layer on the material.

Possible sources of chlorides:

- Exhausted capacity of the ion exchanger or improper maintenance of the system
- If membranes are defective, water from reverse osmosis systems may contain chloride ions
- Use of biocides containing chloride
- Use of drinking water that has been sterilised by „high chlorination“, which may contain chloride ions as well as organic chlorine compounds that are difficult to break down.
- In factories where hydrochloric acid is used: HCl gas may enter from the air
- Input of salts from the ambient air, e.g. salt aerosols near the sea.
- Exposure to gases or vapours containing chlorine or chloride must be avoided; chlorine from disinfectants also leads to pitting corrosion
- Acidic ambient air, water with a high iron content or seawater also lead to corrosion


It is also recommended not to use untreated or unpurified river or cooling tower water that does not comply with the specifications for water content.

The cheapest and simplest solution in case chloride being present is a regular water exchange with thorough cleaning with soft water. (see recommendation above)

All statements made here also apply to bromide ions.

Corrosion resistance of stainless steel at different temperatures and chloride concentrations

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Corrosion

Corrosion is the reaction of a material with chemical substances from its environment („attack agent“). The material is changed, attacked and, in the worst case, even destroyed. Corrosion in aqueous solutions is an electrochemical process.

Corrosion can occur in many different forms, depending on the cause of corrosion. The relevant forms of corrosion are surface corrosion, pitting corrosion, crevice corrosion and aeration corrosion.

Surface corrosion

- Description: Slow and uniform change/removal of the surface.
- Cause: e.g. pH value in the water circuit too low
- Consequences: optical change of the surface, e.g. dulling

Pitting corrosion

Description: Locally increased corrosion without mechanical stress
Cause: Usually increased content of chloride ions in the water circuit, resulting in local destruction of the protective oxide layer on the surface.
Consequences: Pinprick-like or crater-shaped change of the surface.
Hazard potential: Critical, as leaks can occur suddenly.

Crevice corrosion

- Description: Locally increased corrosion in narrow crevices
- Cause: Inhibited diffusion causes a lack of oxygen in narrow crevices.
- In these areas, the pH value can drop sharply and metal ions from the surface are dissolved
- Consequences: Alteration of the material in places that are difficult to access
- Potential danger: possible late detection of corrosion because corroded areas cannot be seen or are difficult to see

Aeration corrosion

- Description: Locally increased corrosion below the water level
- Cause: different concentration of oxygen in deeper layers of the water
- Caused by only partial filling of tanks and lack of circulation
- Consequences: Local corrosive change of the surface

Decalcification

If there are visible limescale deposits, the temperature control unit must be descaled. To do this, empty the temperature control unit and, depending on the degree of calcification, prepare 5 - 10 % citric or acetic acid solution and add to the water circuit. Allow a maximum reaction time of 30 minutes (preferably with circulation) and then rinse thoroughly with soft water (at least 3 times), all areas that have been in contact with the acid.

Dispose of in accordance with company and legal regulations.

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Ambient Temperature Range

is the temperature range which can be attained at an ambient temperature of 20 °C by the circulator alone and with the exclusive use of electrical energy. The operating temperature, that may only be reached by using auxiliary devices, is indicated in brackets. In the case of a heating circulator the working temperature begins above room temperature (as a result of the energy introduced by the pump and the effective insulation) and ends at the upper limit of the operating temperature. The WTR of a cooling circulator begins with the lowest operating temperature of the unit and finishes with the upper temperature at which the refrigeration machine can permanently operate.

Viscosity

indicates the flow characteristics of the heat transfer fluid and it is very dependant on temperature. When operating a unit, the viscosity has a decisive influence on the temperature stability and the pump capacity.

Firepoint

is the temperature at which bath fluids continue to burn after ignition (> 5 sec). According to standard, the maximum working temperature must always be 25 °C below the fire point.

Flashpoint

is the lowest temperature at which emerging vapors can temporarily ignite (< 5 sec). The highest operating temperature specifications for bath fluids are always below the flash point and therefore exclude these risks.

Boiling Point

is the temperature at which a liquid begins to boil.

Ignition Temperature

is the lowest temperature at which the bath fluid ignites spontaneously and continues to burn without heat supply.

Inspired by **temperature** designed for you

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