

# Unistat® 425w

## Cooling a 5-litre HWS reactor

### Requirement

This case study looks at the performance of a Unistat 425w cooling a 5-litre glass reactor from 20 °C to -20 °C and then to  $T_{min}$  under "process" control.

### Method

The Unistat 425w is connected to the 5-litre HWS glass reactor using two insulated metal 1.5-metre hoses. The reactor is filled with 3.75 litre of "M90.055.03", a silicon based HTF.

### Results

The jacket cools quickly to -33 °C to pull the process to its new set-point from 20 °C to -20 °C (40 K) within 15 minutes (ramp rate > 2.6 K/min.).

Once temperatures are stable a set-point of -40 °C is entered. After approximately 22 minutes the jacket temperature asymptotes at -40 °C with a corresponding end-process temperature of -39 °C.

### Setup details

Unistat® 425w & 5-litre HWS reactor

- Temperature range: -40...250 °C
- Cooling power: 2.8 kW @ 250...100 °C  
2.5 kW @ 0 °C  
1.9 kW @ -20 °C  
0.2 kW @ -40 °C
- Heating power: 2.0 kW
- Hoses: 2x1.5 m; M38x1.5 (#6656)
- HTF: DW-Therm (#6479)
- Reactor: 5-litre jacketed glass reactor
- Reactor contents: 3.75 litre M90.055.03 (#6259)
- Reactor stirrer speed: 200 rpm
- Control: process

