



### Setup details

Unistat<sup>®</sup> 815 & Miniplant DDPS 6 L

Temperature range:	-85250 °C
Cooling power:	1.5 kW @ 25020 °C
	1.4 kW @ 40 °C
	1.2 kW @ -60 °C
	0.2 kW @ -80 °C
Heating power:	2.0 kW
Hoses:	2x1 m; M38x1.5 (#9616)
HTF:	M90.055.03
Reactor:	6-litre insulated jacketed
	glass reactor
Reactor content:	4 litre M90.055.03
Stirrer speed:	~ 200 rpm
Control:	process

# Unistat<sup>®</sup> 815

Controlling a 6-litre DDPS jacketed glass reactor

#### Requirement

Chemistry at cryogenic temperatures raises questions on the level of control that can be expected at low temperatures. This case study looks at the ability of a Unistat 815 to cool and control a "DDPS" 6-litre glass reactor to low temperatures.

## Method

The DDPS reactor was connected to the Unistat 815 using two 1-metre long M38x1.5 insulated hoses. The HTF used was Huber's "M90.055.03", a silicon based HTF.

#### Results

The second and third curves demonstrate the pin point and repeatable control made possible through "TAC" control as the Unistat 815 ramps the process temperature between 20 °C and -50 °C and back then from 20 °C to -30 °C.

