



Setup details

Unistat® 510w & Chemglass 50-litre reactor

Temperature range: -50...250 °C

5.3 kW @ 250...0 °C Cooling power:

2.8 kW @ -20 °C 0.9 kW @ -40 °C

6.0 kW Heating power:

Hoses: 2x1.5 m; M38x1.5

(#6659)

HTF: DW-Therm (#6479) 50-litre Chemjacketed Reactor:

> glass reactor (un-insulated)

37 litre M90.055.03 Reactor content:

(#6259)

80 rpm Stirrer speed: Control: process

Unistat® 510w

Cooling a Chemglass 50-litre jacketed glass reactor from 120 °C to -30 °C

Requirement

This case study looks at the speed of response to cool a Chemglass 50-litre jacketed glass reactor to -30 °C from 120 °C (150 K).

Method

The Unistat and reactor were connected using two 1.5 m insulated metal hoses. The reactor was filled with 37 litre of "M90.055.03", a Huber supplied silicon based HTF.

It can be seen that the jacket rapidly cools to -42 °C pulling the process temperature towards its new set-point before heating slightly to guide the process to -30 °C with negligible under-shoot.

