

## Unistat® 910w

Heating a Diehm 100-litre jacketed glass reactor from -80 °C to 20 °C

### Requirement

This case study shows the effectiveness of a Unistat 910w in heating a Diehm 100-litre jacketed glass reactor from -80 °C to 20 °C.

### Method

The Unistat and reactor are connected using two 1.5-metre insulated metal hoses. The reactor is filled with 75 litre of "M90.055.03", a Huber supplied silicon based HTF.

### Results

It can be seen that the jacket temperature ramps rapidly through 165 K (-80 to 85 °C) within 1 hour then ramping quickly down as the process approaches its target to guide the process temperature to its new set-point. Though oversized for a Unistat 910w (designed for efficient operation on reactors to a maximum of 50 litre), the speed and accuracy of the control is evident in the graphic below.

### Setup details

Unistat® 910w & Diehm 100-litre reactor

Temperature range: -90...250 °C  
Cooling power: 5.2 kW from 250 °C to -20 °C

4.7 kW @ -40 °C  
3.1 kW @ -60 °C  
0.9 kW @ -80 °C

Heating power: 6.0 kW  
Hoses: M38x1.5; 1x 2m #6657; 1x1m # 6655, VPC Bypass installed

HTF: M90.055.03 (#6259)

Reactor: 100-litre Diehm un-insulated jacketed glass reactor

Reactor content: 75 litre M90.055.03

Stirrer speed: 410 rpm

Control: process

