

### Unistat® Tango

Unistat Tango controlling a 6l QVF reactor

**Requirement**

The graphics illustrate the performance of Unistat Tango working with a 6l QVF reactor.

**Method**

The Unistat and reactor are connected using two 1,5-metre insulated metal hoses. The reactor is filled with SilOil M40.165/220.10.

**Setup details**

- Temperature range: -45 ... +250 °C
- Cooling power: 0,70 kW @ 0 °C  
0,40 kW @ -20 °C  
0,40 kW @ -30 °C
- Heating power: 1,5 kW
- Hoses: M24 x 1,5 m
- HTF: SilOil M40.165/220.10
- Reactor: 6-litre glas reactor
- Reactor content: 5l M40.165/220.10
- Reactor stirrer speed: 260 rpm
- Control: Process





### Results

**Temperature control of the reactor between 0 °C and +20 °C:**


It can be seen from the graphic how quickly the jacket ramps creating a wide difference in temperature between the jacket and process in the initial cool down phase. Around 33 minutes after the start 0 °C could be reached as process temperature.

In the heat up phase the Unistat Tango takes 21 minutes to heat the 6-litres reactor from 0 °C to +20 °C. The heating rate of 0,95 K/min can be seen on the process temperature curve.



**Temperature control of the reactor in a temperature range from +20 °C to +150 °C:**

The graphic shows the time taken to cool down and heat up the process in a temperature range from +20 °C to +150 °C. The table given below shows the various time taken to cool down and heat up the process in a different temperature ranges.



<u>Ramps / Set point</u>	<u>Time to reach Set-point</u>
+20°C to +50°C	25 min
+50°C to +20°C	32 min
+20°C to +100°C	51 min
+100°C to +20°C	55 min
+20°C to +150°C	1h 22 min
+150°C to +20°C	1h 12 min