



Ministat[®] 240

Ministat[®] 240 cycling a 5-litres glass vacuum insulated reactor

Requirement

This Case Study demonstrates the temperature control capabilities of the process temperature and also the minimum achievable process temperature when a Ministat 240 is connected with a 5-liter Asahi glass vacuum insulated reactor.

Method

The 5-litres Asahi glass vacuum insulated reactor was connected to Ministat[®] 240 using 1-meter metal insulated hoses. The thermofluid used in the system was "M60.115/200.05". "Process" control was carried out via a Pt100 sensor located in the "process" mass. Stirrer speed was set to 150 rpm.


Setup details

Temperature range: Cooling power:	-45°C+200°C 0.60 kW @ +20°C 0.55 kW @ 0°C 0.35 kW @ -20°C
Heating power:	2.0 kW
Hoses:	2*1 m metal insulated
HTF:	M60.115/200.05
Reactor:	Asahi 5-litres glass vacuum insulated
Reactor content:	4 M60.115/200.05
Stirrer speed:	150 rpm
Control:	process
Amb. temperature:	+25°C

Results

1. Performance:

The graphic shows the tight and rapid control as the Ministat 240 ramps to each new set-point. The Ministat[®] 240 needs 57 minutes to cool down the reactor from +20°C to -20°C, then 52 minutes to heat up the reactor from -20°C to +100°C and 64 minutes to cool down the reactor from +100°C to +20°C.





2. Lowest achievable temperature (Tmin):

Once stable at +20°C under the "Process" control, a set point of -40°C is entered. The graphic shows that the lowest temperature achieved in a 1-liter Chemglass jacketed reactor was -29.6°C.

