



Ministat[®] 230

Ministat® 230 cycling a 2000-mL round bottomed reactor

Requirement

This Case Study demonstrates the temperature control capabilities of the process temperature when a Ministat 230 is connected with a 2000mL Chemglass round bottomed reactor.

Method

The 2000-mL Chemglass round bottomed reactor was connected to Ministat® 230 using Tygon A-60-G hoses. The thermofluid used in the system was "M40.165/220.10", reactor was filled with 1000 mL Isopropyl Alcohol. "Process" control was carried out via a Pt100 sensor located in the "process" mass. Magnetic stirrer speed was set to 200 rpm.

Setup details

Temperature range: -40°C...+200°C Cooling power:

	0.38 kW @ 0°C
	0.25 kW @ -20°C
Heating power:	2.0 kW
Hoses:	Tygon A-60-G
HTF:	M40.165/220.10
Reactor:	Chemglass 2000-mL
Reactor content:	round bottomed reactor
Reactor content.	Isopropyl Alcohol
Stirrer speed:	200 rpm
Control:	process

0.42 kW @ +20°C

Results

Performance:

The graphic shows the control of this simple bench-top application.

Start T	End T	Time Taken	Av. Ramp Rate
+20°C	-27.8°C	T-min	XXXX
-27.8°C	+20°C	23 minutes	2.0 K/min
+20°C	+50°C	23 minutes	2.0 K/min
+50°C	+20°C	23 minutes	2.0 K/min

