



Piccolo 280 OLÉ

Rapid heating and cooling of a VMA 0.25 I jacketed stainless steel container

Requirement

This case study demonstrates the heating and cooling performance of a Piccolo 280 OLÉ connected to a VMA 0.2 liter non insulated jacketed stainless steel container.

Method

The Piccolo and container are connected using two 1.0 meter insulated silicone hoses. The container is filled with 0.2 liter of water. Reactor content is stirred by a magnetic stirrer.

Setup details

Unit

Temperature range: +4°C...+70°C Cooling power: 0.28 kW @ +20°C 0.62 kW @ +20°C Heating power:

Container

Manufacturer: VMA Volume: 0.25 liter

Type: Non insulated jacketed stainless steel container

Test conditions

Water 0,2 l water Reactor content:


Stirrer speed: 100 rpm (Magnetic stirrer)

Internal (jacket) Control:

Results

1. Stability


The graphic shows the temperature stability of +/-0.1K at +20°C.





2. Perfomance

Start T	End T	Time Taken	Av. Ramp Rate
+20°C	+60°C	8 Minutes	5.0 K/min



Start T	End T	Time Taken	Av. Ramp Rate
+20°C	+5°C	15 Minutes	1.0 K/min

