



# Piccolo 280 OLÉ

The temperature stability evaluation in case of using a nylon-made fluid pad

#### Requirement

This case study demonstrates the stability of temperature performance of a Piccolo 280 OLÉ connected to a flexible nylon jacket which wraps a glass beaker.

#### Method

The Piccolo and the jacket are connected using two 1.0-meter insulated silicone hoses. The jacket is wrapped around a beaker contains 500ml of water which substitute viscometry sample.

The beaker content is stirred by a magnetic stirrer and the viscometer spindle is immersed in it.

#### Setup details Unit

 Temperature range:
 +4°C...+70°C

 Cooling power:
 0.28 kW @ +

 Heating power:
 0.62 kW @ +

 Hoses:
 2x1.0 m; ID 6

Glass beaker Sample content:

Туре:

Test conditions

HTF:WateStirrer speed:100 rRoom temperature :24°C

0.28 kW @ +20°C 0.62 kW @ +20°C 2x1.0 m; ID 6.0 mm silicone

0,5 l water Glass beaker in a flexible nylon jacket

Water 100 rpm (Magnetic stirrer) : 24°C

## Results

### 1. Stability

Temperature stability of sample was +/-0.4 K at +25°C. This is ideal solution for temperature control of sample in beaker.

