

## Case Study CS 1301

# **Unistat P915w**

Unistat P915w controls a 80 liter De Dietrich reactor

#### Requirement

This case study demonstrates the performance of the Unistat P915w to control the process temperature during simulated exothermic reactions at +20°C and -40°C in "real" ambient conditions. Case study also demonstrates the lowest achievable temperature in the process along with cool down & heat up from +20°C to -60 °C to +20°C.

The tables and the graphics show the responsive, tight and stable control with the jacket temperature being continually adjusted to return and hold the process temperature at the set-points as the thermal load generated by the immersion heater is suddenly changed.

#### Method

To simulate the exothermic reactions, a 600 watt immersion heater was placed inside the reaction mass. The heat output was controlled by a regulator with the results recorded using Huber's "Service software".

#### Setup details

lemperature range.	-90
Heating power:	6.0
Hoses:	2 x
HTF:	MS
Reactor:	De
Reactor content:	60
Stirrer speed:	85
Control:	pro
Amb. temperature:	+2


-90°C...+250°C 6.0 kW 2 x M30 Metal Insulated M90.170.02 De Dietrich 80 liter 60 l DW-Therm 85 rpm process +20°C

### Results

#### 1. Performance:

Controlling and regulating temperature at -40°C with simulated exothermic reactions of 600 Watts (516 kcal/h) and 300 Watts (258 kcal/h).


[	Set Point	Exotherm	ΔT (b/w Process & Jacket)
	-40°C	600W	25K
	-40°C	300W	10K








Cooling from  $+20^{\circ}$ C to  $-60^{\circ}$ C.




Heating from -60°C to +20°C.

Start (°C)	Set Point(°C)	Time (in jacket)	Time (in process)
-60°C	+20°C	50 min	1 h 20 min





Controlling and regulating temperature at +20°C with simulated exothermic reactions of 600 Watts (516 kcal/h) and 300 Watts (258 kcal/h).



2. Lowest achievable temperature in the process:

Start (°C)	T-min (in Process °C)	Time (in process)	T-Min (in Jacket °C )	Time (in process)
+20°C	-88°C	6 h	-90 °C	3 h 30 min

