

Unistat® petite fleur®

Petite Fleur® heating and cooling a 0.5-litre reactor between 20 °C and -30 °C

Requirement

This case study looks at the repeatability of control as the Unistat Petite Fleur cycles the process temperature of a 0.5-litre un-insulated glass reactor.

Method

The Unistat Petite Fleur is connected to the reactor with two insulated metal 1-metre hoses. The Petite Fleur is then programmed to cycle between low and high temperatures.

Results

The Process is cooled rapidly from 20 $^{\circ}\text{C}$ to -30 °C hitting exactly -30 °C with no overshoot. The heating ramp is very fast resulting in a 2 °C overshoot before maintaing perfect stability at 20 °C.

Setup details

Petite Fleur® & Schlee GmbH

Temperature range: -40...200 °C Cooling power: 0.48 kW @ 200...0 °C

0.27 kW @ -20 °C

Heating power: 1.5 kW

Hoses: 2x1 m; M16x1 (#9608)

HTF: Ethanol

Reactor: 0.5-litre un-insulated

glass reactor

Reactor content: 375 ml M90.055.03

(#6259)

Stirrer speed: 160 rpm Control: process



