



## Setup details

Reactor content:

Unistat® 910w & Chemglass reactor

Temperature range: -90...250 °C

5.2 kW @ 250...-20 °C Cooling power:

6.0 kW Heating power:

Hoses: 2x1.5 m; M30x1.5 (#6386) HTF: DW-Therm (#6479) Reactor: 50-litre un-insulated

jacketed glass reactor 35 litre M90.055.03

(#6259)

80 rpm Stirrer speed: Control: process

# Unistat® 910w

Heating and cooling a Chemglass 50-litre un-insulated glass reactor

## Requirement

The graphic shows the performance of a Unistat 910w connected to an un-insulated 50-litre reactor with M30x1.5 hoses.

### Method

The Unistat and reactor are connected using two 1.5-metre insulated metal hoses. The reactor is filled with 35 litre of "M90.055.03", a Huber supplied silicon based HTF.

### Results

From a starting temperature of 20 °C, the purpose of the test is to illustrate the ramp speed of the process temperature in response to setpoint of 100 °C and then back to 20 °C. The process temperature reaches 100 °C from 20 °C within 53 minutes (ramp rate > 1.3 K/ min.). The cooling ramp-rate is > 1.7 K/min. reaching 20 °C from 100 °C in 47 minutes.

