

Unistat® 930w

Cooling a Diehm 100-litre glass reactor to

Requirement

This case study shows the performance of a Unistat 930w connected to a 100-litre Diehm glass reactor cooling from 20 °C to T_{min} under "internal" (jacket) control.

Method

The Unistat and reactor are connected using two 1.5-metre insulated metal hoses. The reactor is filled with 75 litre of "M90.055.03", a Huber supplied silicon based HTF.

Results

The test is run for 2 hours. The initial "internal" (jacket) ramp rate of over 8.5 K/min. rapidly cools the jacket from 20 °C to -70 °C in approximately 12 minutes with a corresponding process ramp rate averaging 2.2 K/min.

After 2 hours the minimum internal temperature reached is -85 °C with a corresponding process temperature of -71 °C though it is continuing to cool towards the "internal" (jacket) temperature.

Setup details

Unistat® 930w & Diehm reactor

Temperature range: -90...200 °C

20 kW @ 0...-40 °C Cooling power:

15 kW @ -60 °C 5 kW @ -80 °C

24 kW Heating power:

Hoses: 2x1 m; M38x1.5 (#6656) HTF: DW-Therm (#6479) 100-litre un-insulated Reactor:

glass reactor VPC Bypass installed

75 litre M90.055.03 Reactor content:

(#6259)

Stirrer speed: 400 rpm Control: internal



