

# Unistat® 910w

## Cooling a Buchi Glas Uster 20-litre reactor to T<sub>min</sub>

### Requirement

This case study is to find the minimum temperature that a Unistat 910w can take the jacket and the resultant process temperature of a Buchi Glas Uster 20-litre glass reactor.

### Method

The Unistat and reactor are connected using two 1.5-metre insulated metal hoses. The reactor is filled with 15 litre of "M90.055.03", a Huber supplied silicon based HTF.

### Results

Initially the jacket temperature ramps at a rate of 11 K/min. and begins to asymptote between -60 °C and -65 °C and finally bottoming out at -82 °C.

The minimum resultant process temperature is around -75 °C. If the test was allowed to run on this would probably have cooled another degree or so.

### Setup details

Unistat® 910w & Buchi Glas Uster reactor

- Temperature range: -90...250 °C
- Cooling power: 5.2 kW @ 250...-20 °C  
4.7 kW @ -40 °C  
3.1 kW @ -60 °C
- Heating power: 6.0 kW
- Hoses: 2x1.5 m; M38x1.5 (#6656)
- HTF: DW-Therm (#6479)
- Reactor: 20-litre un-insulated jacketed glass reactor (#6259)
- Reactor content: 15 litre M90.055.03
- Stirrer speed: 70 rpm
- Control: internal

