

# Unistat Grande Fleur

**Unistat Grande Fleur controls the process temperature in non-insulated 2l glass jacketed reactor from Radleys**



**Requirement**

This case study demonstrates the ability of the Unistat Grande Fleur to control the process temperature in non-insulated 2l glass jacketed reactor from Radleys.

**Method**

The Unistat Grande Fleur was connected to a 2l Radleys non-insulated glass jacketed reactor via 2 x 1m metal insulated tubes. The HTF used was Huber's M40.165/220.10 and the process mass simulated with 1.5l of Huber's M40.165/220.10 silicon oil.

Under "Process Control" from a Pt100 (Teflon covered) located in the process mass, different set points were entered and the performance of the Unistat Grande Fleur was recorded using Huber's service software and recorded onto a USB thumb drive inserted in the USB interface on the Pilot ONE controller.

The agitator speed was set to 250rpm.

**Setup details**


- Temperature range: -40°C...+200°C
- Heating power: 1.8 kW
- Cooling power: 0.6 kW @ +100°C  
0.6 kW @ 0°C  
0.35 kW @ -20°C  
0.2 kW @ -30°C
- Hoses: 2 x M24x1m Metal Insulated
- HTF: M40.165/220.10
- Reactor: Radleys 2l
- Reactor content: 7l M40.165/220.10
- Control: process
- Stirrer speed: 250 rpm
- Amb. temperature: +25°C

**Results**

**1. Temperature Control: from +20°C to +100°C**

This test demonstrates the speed and accuracy that the Unistat Grande Fleur controls the process temperature from +20°C to +100°C and back to +20°C.

Start (°C)	End (°C)	Approximate time (min)	Average Ramp Rate (K/Min)
+20	+100	27	2.96
+100	+20	33	2.42



## 2. Temperature Control: -20°C to +100°C to -20°C

This test demonstrates the speed and accuracy that the Unistat Grande Fleur controls the process temperature from -20°C to +100°C and cooling down to -20°C.

Start (°C)	End (°C)	Approximate time (min)	Average Ramp Rate (K/Min)
-20	+100	36	3.33
+100	-20	66	1.82

