

Unistat Tango

Unistat Tango controls the process temperature in vacuum insulated 10l glass jacketed reactor from Asahi

Requirement

This case study demonstrates the ability of the Unistat Tango to control the process temperature in vacuum insulated 10l glass jacketed reactor from Asahi.

Method

The Unistat Tango was connected to a 10l Asahi vacuum insulated glass reactor via 2 x 1.5m metal insulated tubes. The HTF used was Huber's M40.165/220.10 and the process mass simulated with 6l of Huber's M40.165/220.10 silicon oil.

Under "Process Control" from a Pt100 (Teflon covered) located in the process mass, different set-points were entered and the performance of the Unistat Tango was recorded using Huber's service software and recorded onto a USB thumb drive inserted in the USB interface on the Pilot ONE controller.

Setup details

Temperature range: -45°C...+250°C

Heating power: 3.0 kW

Cooling power: 0.7 kW @ +100°C

0.7 kW @ 0°C

0.4 kW @ -20°C

0.06 kW @ -40°C

Hoses: 2 x M30 x 1.5 m Metal Insulated

HTF: M40.165/220.10


Reactor: Asahi 10l

Reactor content: 6l M40.165/220.10

Control: process

Stirrer speed: 100 rpm


Amb. temperature: +25°C



Results

1. Lowest achievable temperature (Tmin):

The graphic below demonstrates a minimum achievable process temperature of -34.35°C with a corresponding jacket temperature of -35.62°C.



2. Temperature Control

This test demonstrates the speed and accuracy of the Unistat Tango to control the process temperature from +20°C to -30°C, then to +100°C and back to +20°C.

Start (°C)	End (°C)	Approximate time (min)	Average Ramp Rate (K/Min)
20	-30	147	0.34
-30	+100	82	1.58
+100	+20	108	0.74

